




Instructions Slim Light Pro 3000

We want to thank you for your purchase and welcome you to Slim light pro where you can get a great quality machine for your personal use or your small business at an affordable price. Below we will explain how to properly use your device. If you have any questions not answered here feel

free to contact us directly.

Below is a picture of the Slim Light pro 3000 blanket. This was created for those who want simplicity and total body coverage without having to individually wrap each body part. The blanket can be placed on a flat surface and then the client can simply lay on top of it. You can also turn the blanket on and lay it over the top of the person's body. Some also purchase 2 blankets and lay one on top and one below the client. When laying the blanket on a surface make sure it is a firm surface and not too soft. A massage table is ideal to lay the blanket on and a pillow added for comfort.



1. Newest state of the art diodes

Older original diodes needed to increase the level of heat output to put out more brightness and penetration into the body. The problem with this is to achieve deeper penetration you could possibly burn the skin or at least make the session very uncomfortable. The benefits achieved with the lights have nothing to do with the heat output. It's actually quite the opposite. Our new diodes allow us to put out even more power than the older diodes and do this at a cooler temperature. Not only is it more comfortable, it's also more effective for benefits in the body. When the body heats up it will go into a defensive mode and not be as open to absorb as much light to receive all of its many benefits. With your new Slim Light Pro System you are able to achieve maximum benefits during all of your sessions

2. Setting up your system

When setting up your blanket you will have 3 components. Your blanket, the power pack and cord and the control panel .



Mat



Controller and

power supply


First lay out your blanket on a flat surface. Towards the top right when looking at the blanket with lights up you will notice the power Insert wire. This is where you will insert the control handle. Plug in the cord into a standard 110 outlet then into the other side of the control handle. At this point you will notice a red light on the control panel. This is your on/off button. Hold it down for 3 seconds to power up your device. Now that your blanket is on you will need to set the power and the timer. Settings go between 1 and 5 for power output. Once this is established you can set the timer anywhere between 15 and 90 minute sessions. Once set, hit the power button again to start your session.

The Slim Light Pro 3000 has 3 modes to choose from. Most of the time we recommend using the main setting which will have your LED lights on constantly. The unit also has pulse and fast pulse features. These modes can be used to help with more energy levels in the body; they will also penetrate slightly deeper into the body. To switch to the pulse and fast pulse mode while the session is going on press power once for pulse and again for quick pulse. Everyone is different and will have their own preference so we suggest starting with the constant setting and experiment with the extra settings and see how your body feels from using them. Some report more energy and some say they may help more with the deeper penetration for pain in the body. The steady light setting keeps the lights and energy flowing constantly into the body. For overall benefits this is what we recommend for settings.



Specifications

Size	70*27inch
LED Number	4080leds
Wave Length	660nm*850nm
Each Lamp Energy	200mw
Voltage	110V/220V
Material	PVC Cover
Magnetic Field	0 uT
Electric Field	24 V/m
Package Size	89*29*18
Weight	6.7KG
Function	Weight loss/pain relief/anti -aging/skin rejuvenation



4、Cautions

1. During treatment, please do not look directly at the light source.
2. Check whether the temperature sensitivity of the skin at the treatment site of the patient is normal.
3. If there is a wound on the irradiated site, it should be cleaned first. Fresh scar and skin graft should be used with caution.
4. Acute trauma should not be treated locally within 24-48 hours to avoid aggravating swelling, pain and bleeding.
5. During the treatment, patients should be frequently asked and mainly observed. If the patient feels overheated, flustered, dizzy and other reactions, the doctor should be informed immediately.
6. After sessions customers are advised to take a 30-minute walk or exercise and to drink more water to flush out and detoxify the body. Avoid oily foods and sweets. Make dinner a smaller meal and to reduce intake of caffeine and alcohol.

5. Recommended time our possible benefits

Lipo wrap/ therapy wrap

Treatment Time: Weekly

Body sculpting: 30 min 3-5x

Skin rejuvenation: 15min 3-5x

Wound healing: 15 min 3x

Muscle gain: 15 min 3-7x

Recovery: 15 min 3-7x

Pain: 15 min 7x

Warnings

- *Do not stand on the light pads.
- *Make sure not to lift the device by the cord.
- *When unwrapping the client from the pads make sure no stress is applied to the cord coming out of the pads.
- * Do not use the device on open wounds.
- * When using the pads and blankets make sure the client is laying on a firm surface. This is a wired system and by placing it on a firm surface you will be protecting the device to insure your pads will last a long time.

How many sessions are needed?

A normal course of treatment time is between 8 and 12 sessions on average over a 4 to 6 week period of time. Some will see and feel results faster and for some it may take longer. Be consistent and drink plenty of water while doing sessions. We also recommend for best benefits not eating or drinking anything but water and hour before and an hour after treatments.

Will it work for Everyone?

Not everyone is a candidate. Our experience shows that anyone that has any serious medical conditions should check with their doctor before doing treatments. Though the lights are extremely healing for the body, people with pacemakers, active cancer, acute inflammatory lymph issues, liver or kidney disorders, breastfeeding or pregnancy should refrain from use or check with their doctor to see if they feel this would be ok for them to try.

What is red light therapy?

Red light therapy, also known as photobiomodulation, is a therapy designed to use the power of light at certain wavelengths to stimulate the healing power of the body. This FDA approved treatment is safe and effective for many possible health benefits

How does red light therapy work?

Photobiomodulation is directly pro-metabolic as it induces an increase in ATP (energy) production from the mitochondria. Our mitochondria require light to function properly. Red light supplements environmental light with concentrated doses of red and near infrared light which is the specific spectrum the mitochondria uses to produce ATP. By energizing the cells in the body and giving them the energy to do more efficiently what they are meant and designed to do like they did for us when we were younger.

Warranty

We offer 1 year warranty and online support on all products we sell please feel free to call with any questions.

We hope your Slim Light Pro will bring you many years of health for you and your family. Or perhaps to be a great addition to your business. Either way you decide to use your Slim Light Pro, know that we are here to help you in any way we can. Any questions feel free to reach out and

ask one of our staff members. Thank you again for your purchase. We appreciate your business and look forward to hearing from you!

Contact us

Slim Light Pro

Phone: 1-888-547-5296

Website: www.slimlightpro.com