

Instructions Slim Light Pro 4000

We want to thank you for your purchase and welcome you to Slim light pro where you can get a great quality machine for your personal use or your small business at an affordable price. Here are the instructions for your device. If you have any questions not answered here feel free to contact us with any questions.


Below is a picture of the slim light pro 4000 pad system which comes with 7 pads. One large pad which is usually used to lay down on, but some choose to use this section to wrap around the midsection of the body. Besides the large pad your system will come with 6 other pads. 2 small, 2 medium and 2 large. On each pad you will see a Velcro strap attached. Once the pad has been placed, you can use the Velcro strap to secure the pad to the desired area while doing your treatment.

Newest state of the art diodes

The original diodes that were used for LED treatment needed to increase the level of heat output to increase brightness and penetration into the body. The problem with this is to achieve deeper penetration you could possibly burn the skin or at least make the sessions very uncomfortable. The benefits achieved with LED lights have nothing to do with the heat output. It's actually quite the opposite. Our new diodes allow us to put out more than double the power output one each LED at a cooler temperature. Not only is this method more comfortable, it's also more effective for benefits in the body. When the body heats up it will go into a defensive mode and will not be as open to absorb the light and all of its many benefits. Your new Slim Light Pro System is equipped with the latest in LED diode technology ensuring you the greatest results possible.

Specifications


PRODUCT DESCRIPTION



Wavelength	635nm+850nm
Laser Type	Low intensity diode laser
Led Number	4770leds
Laser Pads	7pcs(1 large, 4 medium ,2 small)
Touch Screen Size	8 inch
Package Size	46*46*40
Weight	15KG
Electric Field	24 V/m
Magnetic Field	0 uT
Power	150W
Keywords	Red light system
Material	PVC Cover


Function	Weight loss/pain relief/anti - aging/skin rejuvenation
Voltage	110V/220V

Operation Instruction



Setting up your system

Your Slim Light Pro System will have 7 plugs in the back of the unit. Because the pads are different sizes, different amounts of power are required for each one of your pads. It is very important when setting up the unit that each pad goes in the correct port. Lay out all of the pads on a flat surface and remove the wire ties. Each one of your pads will have a tag on it with a letter S1 S2 M1 etc. all the way up to the large panel or L1. Line up S1 with the S1 hole on the control panel. Notice each plug has a small groove in it when sliding them in. Once it slides into the proper hole, take the silver metal circle on each plug and tighten it on to lock in your plug. Do the same for each plug until all 7 plugs are attached.



Turning on your Slim light Pro

First insert the power cord into the back of your machine.

Turn on the power switch on your unit. Once on, you will see the intro screen light up. Press the next button to go to the control pad interface. Here you will see all the pads, identified by small medium and large. Each pad can be set between 1 thru 5 for power output. You can also control the time from 1 minute up to 60 minutes. Because your unit has full body coverage, you may not want to use all of the pads at the same time. You can set pads you don't want to use to level 0. This will keep those units off during your session. Once your settings are complete, hit the start button on the bottom right to begin your session. A chart with suggested times for different benefits will be listed in this manual. When using this unit, the recommended time for most would be a 20-to-30-minute session. Afterwards to help the body detox we recommend some form of exercise. Many of our clients like to get on a vibration unit for approx. 10 minutes to help activate lymph drainage which works hand in hand with the red light therapy. If you don't have one you can get on a treadmill, walk outside, mini trampoline or do 10 to 50 jumping jacks. Again, the easy way is the vibration unit but other methods are just as effective. Once your session is complete make sure to turn off the main power switch. We recommend wiping down the pads with an alcohol-based cleaning wipe on the plastic part of the pad covers.

Cautions

1. During treatment, please do not look directly at the light source.
2. Check whether the temperature sensitivity of the skin at the treatment site of the patient is normal.
3. If there is a wound on the irradiated site, it should be cleaned first. Fresh scar and skin graft should be used with caution.
4. Acute trauma should not be treated locally within 24-48 hours to avoid aggravating swelling, pain and bleeding.
5. During the treatment, patients should be frequently asked and mainly observed. If the patient feels overheated, flustered, dizzy and other reactions, the doctor should be informed immediately.
6. After treating, the customer is advised to take a 30-minute walk or exercise everyday. And drink more water to flush out and detoxify the body. Urinate. Avoid greasy and oily foods and sweets. Make dinner a smaller meal. Reduce intake of caffeine and alcohol.

Recommended time for possible benefits:

Lipo wrap/ therapy wrap

Treatment Time: Weekly

Body sculpting: 30 min 3-5x

Skin rejuvenation: 15min 3-5x

Wound healing: 15 min 3x

Muscle gain: 15 min 3-7x

Recovery: 15 min 3-7x

Pain: 15 min 7x

Warnings!

*Do not stand on the light pads.

*Make sure not to lift the device by the cord.

*When unwrapping the client from the pads make sure not to apply any stress to the cord coming out of the pads.

* Do not use device on open wounds.

* When using the pads and blankets make sure the client is laying on a firm surface. This is a wired system and by placing it on a firm surface you will be protecting the device to ensure your pads stay protected.

How many sessions are needed?

This will depend on a few factors such as: how many inches you want to lose. Body fat % eating habits and exercise. A normal course of treatment time is between 8 and 12 sessions on average over a 4 to 6 week period of time. Some will see and feel results faster and for some it may take longer. Be consistent and drink plenty of water while doing sessions. We also recommend for best benefits not eating or drinking anything but water and hour before and an hour after treatments.

Will it work for Everyone?

Not everyone is a candidate. Our experience shows that anyone that has any serious medical conditions should check with their doctor before doing treatments. Though the lights are extremely healing for the body, people with a pacemaker active cancer, acute inflammatory lymph issues, liver or kidney disorders, breastfeeding or pregnant should refrain from use or check with their doctor to see if they feel this would be ok for them to try..

What is red light therapy?

Red light therapy, also known as photobiomodulation is a therapy designed to use the power of light at certain wavelengths to stimulate the healing power of the body. This FDA approved treatment is safe and effective for many possible health benefits

How does red light therapy work?

Photobiomodulation is directly pro-metabolic as it induces an increase in ATP (energy) production from the mitochondria. Our mitochondria require light to function properly. Red light supplements environmental light with concentrated doses of red and near infrared light which is the specific spectrum the mitochondria uses to produce ATP. By energizing the cells in the body and giving them the energy, they are able to do what they were meant and designed to do. They are able to behave more like they did for us when we were much younger.

Warranty

We offer 1 year warranty and online support on all products we sell. Please feel free to call with any questions.

We hope your slim light pro will bring you many years of happiness and health for you and your family. Or perhaps be a great addition to your business. Either way you decide to use your device, know that we are here to help you in any way we can. Feel free to reach out and ask one of our staff members any questions you may have. Thank you again for your purchase. We appreciate your business and look forward to hearing from you!

Contact us

Slim light pro

Phone 1-888-547-5296

Website: <https://www.slimlightpro.com>